

BP-200 plus

Automated blood pressure monitor for cardiac exercise testing



SCHILLER
The Art of Diagnostics

BP-200 plus

- ∴ Seamless integration: automated communication with your SCHILLER exercise ECG device reduces the risk of transmission errors.
- ∴ Non-exercise BP mode: allows BP measurements during patient setup and recovery without ECG signal.
- ∴ Verified measurements: Korotkoff sounds are displayed.



BP AND SPO₂ MEASUREMENTS FOR THE EXERCISE TEST ENVIRONMENT

The BP-200 plus blood pressure monitor is designed to give you the confidence to know that mechanical vibration, patient movement and observer variability does not interfere with an accurate reading. Your focus can be where it matters most – on your patient. The hands-free interface and SpO₂ measurement option make the BP-200 plus a comprehensive device for diagnostic test centres. The BP-200 plus is also ideal in research environments because raw data can be saved on an SD card.

SCHILLER K-SOUND ANALYSIS FOR ACCURATE MEASUREMENTS

The BP-200 plus excels where other BP monitors fail, by using SCHILLER K-Sound Analysis (SKA). SKA was developed as a result of SCHILLER's pioneering work in 24-hour ambulatory blood pressure monitoring. It was this experience that provided SCHILLER's research and development team with invaluable insights into the relationship between many of the variables involved in determining the blood pressure of a patient in motion. These variables include, but are not limited to, K-sounds, overall body shape, foot strike characteristics, arm movements, along with treadmill and ergometer noise. With SCHILLER K-Sound Analysis, Korotkoff sounds are distinguished from artefact noise so that blood pressures are measured accurately with greatly reduced false measurements.

AUTOMATED COMMUNICATION WITH SCHILLER EXERCISE ECG SYSTEM

Recording an accurate and reliable BP measurement or SpO₂ value at the correct intervals during an exercise ECG test can be difficult and stressful for clinicians and lab technicians. That is why the BP-200 plus adopts the programmed intervals from the SCHILLER exercise ECG device and therefore makes sure that measurements are always taken at the same point in time during each stage, increasing measurement reliability.

Systolic BP, diastolic BP, mean arterial pressure, SpO₂ (option), peripheral pulse and double product (sys × HR)

5 buttons for easy programming and operation



Display for Korotkoff sounds or trigger waves

SCHILLER BP cuff is designed for best comfort during exercise ECGs





Americas
SCHILLER Americas Inc.
Doral, Florida 33172

North America:
Phone +1 786 845 06 20
Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com

Latin America & Caribbean
Phone +1 305 591 11 21
Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com



Asia
SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 70 90
Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
Alfred Schiller (Beijing) Medical Equipment
Co.,Ltd.
100015 Beijing, China
Phone +86 010 52007020
Fax +86 010 52007020-8016
info@schillerchina.com
www.schillermedical.cn



Croatia
SCHILLER d.o.o.
10000 Zagreb
Phone +385 1 309 66 59
Fax +385 1 309 66 60
info@schillerzg.hr
www.schiller.hr



France
SCHILLER Médical
F-67160 Wissembourg
Phone +33 3 88 63 36 00
Fax +33 3 88 63 36 49
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
F-77608 Bussy St Georges
Phone +33 1 64 66 50 00
Fax +33 1 64 66 50 10
contact@schillerfrance.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
D-85622 Feldkirchen b. München
Phone +49 89 62 99 81 0
Fax +49 89 609 50 90
info@schillermed.de
www.schillermed.de



India
SCHILLER Healthcare India Pvt. Ltd.
Mumbai - 400 059, India
Phone +91 22 6152 3333 / 2920 9141
Fax +91 22 2920 9142
sales@schillerindia.com
www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 843 20 89 / +48 22 647 35 90
Fax +48 22 843 20 89
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
AO SCHILLER.RU
125040 Moscow, Russia
Phone +7 (495) 970 11 33
Fax +7 (495) 956 29 10
mail@schiller.ru
www.schiller.ru



Serbia
SCHILLER d.o.o.
11010 Beograd
Phone +381 11 39 79 508
Fax +381 11 39 79 518
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
2310 Slovenska Bistrica
Phone +386 2 843 00 56
Fax +386 2 843 00 57
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
ES-28232 Las Rozas/Madrid
Phone +34 91 713 01 76
Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es



Switzerland
SCHILLER-Reomed AG
CH-8912 Obfelden
Phone +41 44 744 30 00
Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKIYE
Okmeydani-Sisli - Istanbul
Phone +90 212 210 8681 (pbx)
Fax +90 212 210 8684
info@schiller.com.tr
www.schiller-turkiye.com



United Kingdom
SCHILLER UK
Bellshill, ML4 3PR
Phone +44 1698 744 505
Fax +44 1698 744 474
sales@schilleruk.com
www.schilleruk.com

Device availability in your market is subject to regulatory approval.

